




FAST FACE SHIELDS

USING TRANSPARENCY FILM
OR SHEET PROTECTORS


Surgical Mask full face shield



length of cut = height of mask






angling the slits may be helpful.





Feed mask loops through slits.
Mask worn **INSIDE** shield.

Surgical Mask eye shield only




Feed mask loops through slits.
Center section tucks into mask.

N95 Mask full face shield




angling the slits may be helpful.




Tuck side flaps into straps.
Mask worn **INSIDE** shield.

Glasses/Goggles full face shield




x2 ties - anything for fastening:

- zip ties
- string
- wire twist ties
- elastic bands
- tape
- paper clips



Holes are farther apart
than width of glasses.



Slide arms of glasses through
outermost slits and secure with ties.
Glasses worn **INSIDE** shield.