

FAST FACE SHIELD

MATERIALS

x1 pair of glasses


- Wearer's prescription frames
- safety goggles
- sunglasses with lenses popped out
- readers with lenses popped out

x2 ties

- zip ties
- wire twist ties
- string

x1 sheet of clear plastic

- transparency film
- page protector
- cut 2-liter soda bottle






hole punch OR knife OR scissors

INSTRUCTIONS

1. CUT HOLES in plastic sheet.

Two slits or holes on each side, towards the middle-top



2. FEED TIES through holes

3. ATTACH to sides of glasses frame

* Attach towards the front of the face, allowing the plastic to curve forward and provide a gap of space between the shield and the nose/face.

* For comfort, tie on OUTSIDE of shield, not inside where it could poke+scratch the wearer's skin.

